

Insulin Infusers

Product Name (Manufacturer)	Description/Use
I-Port Injection Port (Patton Medical Devices)	Medication delivery channel: An insertion needle guides a soft cannula under the skin. Once applied, the insertion needle is removed and only the soft cannula remains below the skin, acting as a direct channel into the skin. Device can accommodate 75 injections, and be worn for up to 72 hours.
Insufion Indwelling Cannula (IntraPump Infusion Systems)	Small plastic tube used as delivery channel. Resembles an IV tube but goes under the skin instead of into a vein. A tiny needle is used to insert the cannula, which is held in place with a dressing for 3 to 7 days.