

# Glycohemoglobin Tests

Name (Manufacturer/Distributor)	Comments
<b>AccuBase A1c Test Kit</b> (Diabetes Technologies)	Fingerstick A1C mail-in test kit. FDA approved for over-the-counter use. Sample mailed to CLIA-licensed, CAP reference lab. Lab procedure screens for abnormal hemoglobins such as S, C, and F, checks for abnormal peaks and/or red blood cell disturbances (anemia), and reports an interference-free A1C answer. Positive patient ID, HIPPA compliant (confidential patient reporting). Approved for diabetes screening. Results reported by mail, fax, e-mail (PDF), or via Web site.
<b>A1c At-Home</b> (FlexSite Diagnostics)	Home sample collection kit for mail-in laboratory A1C testing. Dried blood sample from fingerstick is mailed to FlexSite Diagnostics' CLIA-licensed laboratory. Results reported by mail, fax, or electronically. FDA cleared for over-the-counter sales. Test method certified by the National Glycohemoglobin Standardization Program.