

## Combination Blood Glucose/Blood Pressure Monitors

Name (Manufacturer)	Size (inches)	Weight (ounces)	Test Strip Used	Range (mg/dl)	Test Time	Battery	Warranty	How Calibrated	Control Solution	Features
<b>Advocate DUO Speaking Glucose + Blood Pressure Monitor 2 in 1 System 3223A</b> (Diabetic Supply of Suncoast)	2.5 X 3.0 X 4.0 (with wrist cuff)	5.5 with batteries	Advocate	20-600 glucose 0-300 mmHg	7 sec. glucose, approx. 45-50 sec. blood pressure	(2) AAA	1 year	Built-in button	Yes	Test glucose and blood pressure with the same product. Uses 0.7-microliter sample to test. Has a memory of 450 past glucose tests. Calculates 7-, 14-, 21-, 28-, 60- and 90-day averages for past blood glucose results. Gives heart rate. Alternate site testing. Audible results in English and Spanish. Downloadable software and training videos available.
<b>DuoCare</b> (GenExel Sein)	2.25 X 3.25 X 4.75	14	DuoCare	20-600 glucose 40-200 bpm pulse 20-330 mmHg	5 sec. glucose, approx. 30-45 sec. blood pressure	(2) AAA	3 years	Automatic	No	Combined blood glucose and blood pressure monitor. Requires a 0.5-microliter sample size for testing. Stores 150 past blood glucose values in its memory. Wrist monitor stores 60 blood pressure measurements in its memory. Averages 14- and 30-day test results. Monitor comes with lancing device and replacement lancets, test strips, carrying case, logbook.
<b>FORA D10 2-in-1 Wrist Type Talking Blood Glucose + Blood Pressure Monitoring System</b> (Fora Care)	2.48 X 2.95 X 1.57 (without cuff)	5.6 with batteries	FORA D10	20-600 glucose 0-300 mmHg	7 sec. glucose, approx. 40 sec. blood pressure	(2) AAA	2 years	Built-in button	Yes	Tests blood glucose and blood pressure with one device. Requires a 0.7-microliter blood sample to test. Stores 450 past tests in its memory. Ability to average past 7 and 14 days of glucose results. Alternate site testing. Talking function guides test procedure in both English and Spanish. Managed by health care software downloaded to PC or Web.
<b>FORA D15 2-in-1 Arm Type Blood Glucose + Blood Pressure Monitoring System</b> (Fora Care)	5.39 X 3.54 X 2.13 (without cuff)	16.7 with batteries	FORA D15	20-600 glucose 0-300 mmHg	7 sec. glucose, approx. 40 sec. blood pressure	(4) AAA	2 years	Insert code strip	Yes	Tests blood glucose and blood pressure with one device. Uses 0.7-microliter blood sample to test. Stores 352 past tests in its memory. Averages past 7 or 14 days of glucose results. Alternate site testing. Downloadable results to PC or Web.