

# Insulin Pump Supplies

Product Name (Manufacturer/Distributor)	Type of Product	Features
<b>Accu-Chek Rapid-D</b> (Roche Diabetes Care)	Infusion set	90° insertion angle. Surgical stainless steel, 28-gauge needle. 6-, 8-, and 10-mm needle length with 24", 31", or 43" tubing. Simple connect/disconnect. Manual insertion. Various insertion site options. Accommodates people with low vision. Compatible with all standard luer-lock insulin pumps. Hypoallergenic self-adhesive. No insertion device needed.
<b>Accu-Chek Tender</b> (Roche Diabetes Care)	Infusion set	20° to 45° insertion angle. 27-gauge needle. Soft Teflon cannula. 13- and 17-mm cannula lengths with 24", 31", or 43" tubing. Connect/disconnect at site. Kink-resistant tubing. Compatible with all standard luer-lock pumps. Manual insertion. No insertion device needed.
<b>Accu-Chek Ultraflex</b> (Roche Diabetes Care)	Infusion set	90° insertion angle. 27-gauge needle. Teflon cannula with 8- and 10-mm cannula lengths and 24", 31", and 43" tubing. Tapered, rounded edge to reduce risk of snagging and dislodging. Kink-resistant tubing. Simple connect/disconnect at site. Compatible with all standard luer-lock insulin pumps. Accu-Chek Link Assist insertion device now available for Accu-Chek Ultraflex infusion sets for automated insertion.
<b>Cleo 90</b> (Smiths Medical)	Infusion set	90° insertion angle. 28-gauge needle. Soft cannula in 6- and 9-mm lengths and 24", 31", and 43" tubing. Insertion needle automatically retracts into the insertion device. Compatible with all standard luer-lock pumps. Transparent skin adhesive. Insertion device is included as part of set.
<b>Fifty50 Medical Infusion Set</b> (Fifty50 Medical)	Infusion set	90° insertion angle with quick disconnect at site. Flexible, non-metal cannula available in 6- and 9-mm lengths. 23" and 43" tubing. 27-gauge needle. Compatible with all standard luer-lock pumps. Feature to protect site from accidental tubing pulls.
<b>Fifty50 Medical Reservoirs</b> (Fifty50 Medical)	Reservoirs	Model 3050: disposable, single-use insulin reservoir for use with Medtronic MiniMed 504, 506, 507c, and 508 insulin pumps. Model 1850 and 3070: single-use insulin reservoir for use with Medtronic MiniMed Paradigm pump. 1.8- and 3-ml capacity. Design allows Paradigm users to connect any standard infusion set to their pump.
<b>Orbit 90</b> (ICU Medical)	Infusion set	90° insertion angle with 6- and 9-mm cannula lengths and 24", 31", and 41" tubing. 27-gauge needle. Smooth dome shape to reduce snagging. One-handed activation.
<b>Quick Set with Quick-serter</b> (Medtronic Diabetes)	Infusion set	90° insertion angle with 6-mm or 9-mm soft, tapered cannula and 23" and 43" tubing. 27-gauge needle. Compatible with all standard luer-lock pumps. Built-in adhesive. At-site disconnect. Can be inserted manually or via the Quick-serter insertion device.
<b>Silhouette with Sil-serter</b> (Medtronic Diabetes)	Infusion set	Variable insertion angle with 13- and 17-mm soft, tapered cannula. Available with 23" and 43" tubing. 27-gauge needle. Compatible with all standard luer-lock pumps. Built-in adhesive. At-site disconnect. Can be inserted manually or via the Sil-serter insertion device.
<b>Sure-T</b> (Medtronic Diabetes)	Infusion set	90° insertion angle with 6-mm steel needle. 23" tubing that disconnects 4 inches from insertion site for easy detachment from pump. Built-in adhesive. Compatible with Medtronic MiniMed Paradigm insulin pump.